

EZee Life™

Wide Seat Rollator

Model: CH3040



EZee Life™ Products

34 Futurity Gate, Unit 15
Concord, Ontario L4K 1S6



Specifications:

Weight capacity:	286 lbs	(130 kgs)
User Height*:	5' 4" to 6'	(1.6 to 1.83 m)
Handle Height:	30 3/4" to 35"	(78 – 89 cm)
Seat Height:	21 3/4"	(55 cm)
Seat Width:	22"	(56 cm)
Seat Depth:	9.8"	(25 cm)
Wheel Size:	8"	(20 cm)
Overall Width:	23 3/4"	(60.3 cm)
Weight:	20.9 lbs	(9.5 kg)


*User height is approximate and will vary based on user's body proportions

Brake Adjustment:

Over time the brake cables will stretch and require adjustment.

With the brakes NOT applied hold the adjusting barrel (B) in place and rotate the adjusting nut (N) clockwise to increase brake force, rotate counter clockwise.

IMPORTANT: before using the walker again, ensure the brakes are working properly. If the wheels turn while the brakes are applied do not use the walker.



Safety Precautions

1. Consult your doctor, therapist, or other qualified professional to determine if this mobility device meets your specific needs.
2. Before using, make sure all parts are secure and the frame is fully opened.
3. Ensure brakes are engaged before sitting on the rollator seat or standing from the seat.
4. Do not use as transport device to transport people sitting on seat.

Components:



- | | |
|---------------------------|----------------------|
| A. Brake | G. Handgrip |
| B. Push Handle | H. Backrest |
| C. Push Handle Hand Screw | I. Brake Cable |
| D. Backrest Bracket | J. Seat |
| E. Storage Bracket | K. Folding Mechanism |
| F. Rear Wheel | L. Front Wheels |

Assembly Instructions:

A) Unpack Rollator

Remove rollator from its carton and unwrap all the parts


B) Install Push Handles

- 1) Stand the rollator upright and open the frame by pushing down on the folding mechanisms until they are flat.
- 2) Insert the handlebars. Make sure the brake cables are on the outside of the frame
- 3) For proper height adjustment, stand upright with shoulders relaxed and arms at your side. Handles should be at the height where you would wear a watch.
- 4) Secure the handlebars by inserting the hand screw into the hole in the frame. The hand screw should be flush with the frame. Hand tighten only.



C) Install Backrest:

Insert the backrest into the front backrest brackets (1) until the push button click (2) Brake cable should be outside the frame



D) Folding/Unfolding & Installing Basket

Lock the brakes by pushing down on the brake handles until you feel the brakes click into locked position (1)

To fold remove the basket and pull up on the strap (2)

To unfold push down on the folding mechanism until they are fully extended (3) Install the basket and release brakes (4 & 5)



E) Using the Brakes:

Push down on both brake handles to lock the brakes. (1)

Squeeze both brake handles to release the brake lock (2)

CAUTION The rear wheels will lock but it may still slide if the floor surface is slippery.

REMEMBER to push down (on the brake handles) before you sit down!

CAUTION never release the brakes until you are standing stable and in the ready walk position. Do not release the brakes if you are sitting down.

REMEMBER if you are walking and wish to slow down, squeeze the brake handles slightly to put brake pressure on the wheels.

Warranty:

Ezee Life™ products carry a one year warranty against manufacturing defects, faulty materials and workmanship.

In the event of a defect covered by this warranty, we will, at our option, repair or replace the item. Please contact your supplier/dealer for assistance should a defect become apparent.

The warranty does not extend to consequential costs resulting from the fault or defect in the product, in particular freight and travel costs, loss of earnings, expenses or other expenses which may occur.

The warranty does not cover normal wear and tear resulting from normal use or damage caused by accident, misuse, negligence or improper installation or assembly.

The vendor/manufacturer will not accept responsibility for any damages or injury caused by misuse, negligence or non-observance of the instructions supplied.

Imported & Distributed by:

Factory Direct Medical
34 Futurity Gate, Unit # 15
Concord Ontario L4K 1S6
Local: 416-739-8393
Toll Free: 855-235-2400

