

TILT & RECLINE REHAB SHOWER COMMUNE



MODEL: 170

Traditionally, shower/commodes were predominantly constructed from steel or stainless steel, each with its inherent drawbacks—steel being susceptible to rust and stainless steel being both expensive and heavy. In contrast, Ezee Life™ commodes are meticulously fashioned from high-quality aluminum, offering rust resistance, a lighter weight compared to steel, and a more cost-effective alternative to stainless steel.

These Aluminum Shower/commode chairs facilitate seamless client movement, enabling care providers to handle heavier individuals effortlessly. The distinctive design and incorporation of tilt/recline prioritize safety, comfort, and optimal positioning. The ergonomic headrest is adjustable to accommodate the diverse needs of any client. Whether utilizing a sit-to-stand lift, sling lift, or client pivot transfer, this commode/shower chair seamlessly accommodates transfers. The majority of mechanical lifts easily fit under the frame, ensuring secure transfers that can be executed forward-facing or by lifting the arms for side transfers.

Ezee Life commodes boast the softest seat, backrest, and armrests in the market. This not only guarantees client comfort but also provides sustained support during extended periods of sitting. The tactile texture is gentle to the touch, delivering a soothing experience while maintaining a plush and comfortable seat.



SPECIFICATIONS	
Frame Material	Aluminum
Weight	46.5 lbs (21.1 kgs)
Tilt Range	Up to 13°
Recline Range	90° to 111°
Upper Back Angle	15°
Calf Support Range	90° to 50°
Casters	5" Swivel Locking Casters
Seat Width	18" (45.7 cm)
Seat Depth	17" (43.2 cm)
Seat Height Adjustable	21.8" - 23.8" (55.4 - 60.5 cm)
Foot Rest Length	14" to 20" (35.6 - 50.8 cm)
Foot Rest Hanger Angle	60°
Back Height	16" (40.6 cm)
Arm Height	8" (20.3 cm)
Arm Pad Length	11.75" (29.8 cm)
Overall Length	Up to 50.3" (127.762 cm) Dependant on adjustments of the headrest and footrest)
Overall Width	22.2" 56.388 cm)
Weight Capacity	286 lbs (130 kgs)

FEATURES

- Cushioned molded seat, back, & armrest pads
- Commode hole insert
- Flip back padded arms
- Locking swivel casters
- Facilitate sit-to-stand transfers & mechanical transfers with ease
- Adjustable height swing away footrests
- Soft to the touch composite flip up foot plates
- Adjustable seat belt
- Strong tubular aluminum frame
- Combined tilt and recline action



34 Futurity Gate, Unit # 15 Concord Ontario L4K 1S6

wk: 416-739-8393 | Toll Free: 1-855-235-2400

email: Info@factorydirectmedical.com | www.factorydirectmedical.com