

A portable, personal lift that also doubles as a better alternative to a rollator walker, the device has an elevating seat that can descend to floor level, raise to a height of 24" and stop anywhere in between.



- Lifts a fallen person from the floor; assists with daily activities (transfers, bending, toileting, standing)
- Compact, portable form factor for home, clinic, and community settings
- Reduces caregiver strain (and compensation claims)



SKU	RAY-BK	RAY-RD	RAY-BL
Color	Black	Red	Blue
GTIN	00628774072001	00628774072018	00628774072025
MSRP CAD	\$5,733		
Note: MSRP assumes \$200 CAD shipping			

Weight Capacity for User	136 kg / 300 lbs
Height Maximum for User	203 cm / 6' 8"
Width (Unfolded)	61.6 cm / 24.25"
Length	71 cm / 28"
Height	88.9 cm / 35"
Width (Folded)	27.9 cm / 12.5"
Seat Height Range	Floor Level – 61 cm / 24"
Seat Width	48.3 cm / 19"
Handle Adjustment Range (above seat)	27.9 cm / 11.5" – 38.1 cm / 15.5" (up to 100 cm / 39.5" from floor)
Device Weight	16.3 kg / 36 lbs

Battery	• Inspired Energy PH3059QE29 • 30+ Lifts at 300 lbs on a single charge
Charging	• Barrel connector on side cover of device • ~3 hr charge time
Safety Lock / Electric Parking Brake	• Engage when seat is moved up or down • Remain engaged until 'unlock' button is pressed • Can also engage by pressing 'lock' button
Pre-set Seat Height	• Seat memory for user's ideal sit-to-stand height
Rotating Seat	• User can stand or walk inside the frame
Adjustable & Rotating Handles	• Hand adjustable handle height • Rotate handles out of way to enable transfers
Wheel Fenders	• Hand hold while moving onto seat • Can be used as footrests while seated on device