

LIFT CHAIRS

EZEE LIFE PRODUCT DESCRIPTION, FEATURES & SPECIFICATIONS





OUR MISSION

Since 1993, FDM has dedicated itself to developing and providing health care solutions that allow individuals, with physical challenges, the opportunity to safely enjoy the freedom of increased mobility and independence.

OUR VALUES

We continually work closely with industry professionals to stay informed of market trends and our clients changing needs with an ongoing commitment to our four core values:

Continuously Improve—the quality of life for clients through increased comfort, function and mobility.

Be Accountable—for the solutions we provide.

Exceed Expectations—with results that match the needs of the individual.

Be Progressive—by developing, and introducing new products that reflect our level of dedication to improved mobility and independence.

OUR GOAL

With compassion, commitment, and core values intact, we strive to remain a reliable and responsible source of better living solutions.

C O N T E N T S

Features	4
Mars	5
Earth	6
Mercury	7
Saturn	8
Pluto	9
Venus	10
Jupiter	11
Jupiter (Single Motor)	12

LIFT CHAIRS FEATURES



The spring assisted footrest on all models will flip up if it comes in contact with an obstruction when lowering to minimize injury or damage.



Supplied locking swivel casters can be installed on all models, except Earth and Mars which have 3" casters. Ideal for users who may want to move their chair occasionally or rental units.



Our Saturn and Mercury models are upgraded to include a cup holder and adjustable tray for the convenience of the user. Mercury also includes leather upholstery and heat & massage.



The Earth and Mars models are our institutional models that include 3" locking casters, push bars, vertical lift, removable arms and shoulder wings and incontinent proof black upholstery.

FEATURES:

- Lounger style that fully supports the legs (no gap between footrest & seat)
- Removable back for easy handling
- Easy to use hand control
- UL / CUL Approved
- Maximum seat elevation with superb stability in full lift position
- Wide construction provides greater side to side stability
- Top steel base for durability & cross supports
- Engineered bushings & washers at pivot points provide quiet smooth operation & durability
- Direct drive activation synchronizes right & left sides for easy opening
- Low voltage motor, 24 volt
- Battery Back Up



Jupiter

Mars

Saturn

Pluto

Venus

Earth

Mercury

MARS CH4004M (EARTH GREY)

The Mars lift chair is a three position lift chair that, in addition to lifting in a conventional tilt forward manner, will lift vertically to allow side transfers to or from a bed. It comes standard with earth grey incontinent proof upholstery and features four 3" casters for easy repositioning and transportation of users. With detachable armrests and shoulder wings this chair makes side transfers a possibility for those who are unable to stand.

With a Wide/Tall size rating, this chair will accommodate people of wider widths and above average heights. The Mars Lift chair also features a vertical lift function where the entire chair will lift vertically 6" and the arms and shoulder wings can be removed for users who must, or prefer, a side transfer to and from the chair.

- Three position recline
- Four locking 3" swivel casters
- Detachable armrests
- Incontinence proof upholstery
- Push Handle for user transport

MODEL NAMES	MARS
Model numbers:	CH4004M (Earth Grey)
Weight Capacity:	352 lbs (160 kg)
Seat Width:	21.1" (53.6 cm)
Seat Depth:	21.1" (53.6 cm)
Back Height:	26.9" (68.3 cm)
Seat Height:	19.9" (50.5 cm)
Overall Width:	30" (76.2 cm)
Overall Depth:	68.3" (173.5 cm)
Overall Height (down):	41.7" (105.9 cm)
Overall Height (raised):	57.3" (145.5 cm)
Max. Vertical Lift:	6" (15.5 cm)
Box 1	88CM x 97CM x 72CM
Box 2	82CM x 83CM x 38CM
Net weight	60 kg (box1) 6.5 kg (box 2)
Gross Weight	64 kg (box1) 8.5 kg (box 2)



LIFT CHAIRS (2 MOTORS)

EARTH CH4012M (EARTH GREY)

The Earth lift chair is a three position lift chair that, in addition to lifting in a conventional tilt forward manner, will lift vertically to allow side transfers to or from a bed. It comes standard with earth grey incontinent proof upholstery and features four 3" casters for easy repositioning and transportation of users. With detachable armrests and shoulder wings this chair makes side transfers a possibility for those who are unable to stand.

With a Wide/Tall size rating, this chair will accommodate people of wider widths and above average heights. The Earth lift chair also features a vertical lift function where the entire chair will lift vertically 6" and the arms and shoulder wings can be removed for users who must, or prefer, a side transfer to and from the chair.

- Three position recline
- Four Locking 3" Casters
- Detachable armrests
- Push Handle for user transport
- Incontinence proof upholstery

MODEL NAMES	EARTH
Model numbers:	CH4012M (Earth Grey)
Weight Capacity:	352 lbs (160 kg)
Seat Width:	20.7" (52.6 cm)
Seat Depth:	19.9" (50.5 cm)
Back Height:	28.5" (72.4 cm)
Seat Height:	21.1" (53.6 cm)
Overall Width:	29.3" (74.4 cm)
Overall Depth:	68.6" (174.2 cm)
Overall Height (down):	45.2" (114.8 cm)
Overall Height (raised):	58.5" (148.6 cm)
Max. Vertical Lift:	6" (15.2 cm)
Box 1	88CM x 97CM x 72CM
Box 2	79CM x 73CM x 39CM
Net weight	63 kg (box1) 9.8 kg (box 2)
Gross Weight	67 kg (box1) 11.8 kg (box 2)



MERCURY CH4011M (DARK BROWN) CH4013 M (TAN)

The Mercury lift chair is an infinite position lift chair available with dark brown or tan genuine leather upholstery. While it comes fitted with rubber pads to prevent floor damage and provide stable positioning, the 4 supplied locking casters can be installed to ease repositioning the chair. Additionally, the Mercury comes equipped with an adjustable tray and a cup holder for user convenience.

With a Medium/Short size rating, this chair will accommodate people of average widths and shorter heights.

- Infinite position recline
- Four locking swivel casters included
- Genuine leather upholstery
- Manual tray and cup holder
- Heat and Massage

MODEL NAMES	MERCURY
Model numbers:	CH4011M (Dark Brown) CH4013M (Tan)
Weight Capacity:	352 lbs (160 kg)
Seat Width:	19.5" (49.5 cm)
Seat Depth:	19.1" (48.5 cm)
Back Height:	26.9" (68.3 cm)
Seat Height:	19.1" (48.5 cm)
Overall Width:	28.9" (73.4 cm)
Overall Depth:	69.8" (177.3 cm)
Overall Height (down):	42.5" (108 cm)
Overall Height (raised):	55" (139.7 cm)
Max. Vertical Lift:	n/a
Box 1	81CM x 77CM x 65CM
Box 2	76CM x 67CM x 33CM
Net weight	45 kg (box1) 6 kg (box 2)
Gross Weight	47 kg (box1) 8 kg (box 2)



LIFT CHAIRS (2 MOTORS)

SATURN CH4005M (BLUE) CH4006M (BROWN)

The Saturn lift chair is an infinite position lift chair available with brown or blue upholstery. While it comes fitted with rubber pads to prevent floor damage and provide stable positioning, the 4 supplied locking casters can be installed to ease repositioning the chair. Additionally, the Saturn comes equipped with an adjustable tray and a cup holder for user convenience.

With a Wide/Average size rating, this chair will accommodate people of wider widths and average height.

- Infinite position recline
- Manual tray and cup holder
- Four locking swivel casters included

MODEL NAMES	SATURN
Model numbers:	CH4005M (Blue) CH4006M (Brown)
Weight Capacity:	352 lbs (160 kg)
Seat Width:	21.1" (53.6 cm)
Seat Depth:	19.5" (49.5 cm)
Back Height:	28.5" (72.4 cm)
Seat Height:	19.5" (49.5 cm)
Overall Width:	31.2" (79.2 cm)
Overall Depth:	73.3" (186.2 cm)
Overall Height (down):	43.3" (110 cm)
Overall Height (raised):	58.5" (148.6 cm)
Max. Vertical Lift:	n/a
Box 1	81CM x 83CM x 79CM
Box 2	n/a
Net weight	54 kg
Gross Weight	60 kg



PLUTO CH4007M (BLUE) CH4008M (BROWN)

The Pluto lift chair is an infinite position lift chair available with brown or blue upholstery. While it comes fitted with rubber pads to prevent floor damage and provide stable positioning, the 4 supplied locking casters can be installed to ease repositioning the chair.

With a Medium/Average size rating, this chair will accommodate people of medium width and average height.

- Infinite position recline

- Four locking swivel casters included

MODEL NAMES	PLUTO
Model numbers:	CH4007M (Blue) CH4008M (Brown)
Weight Capacity:	352 lbs (160 kg)
Seat Width:	19.5" (49.5 cm)
Seat Depth:	19.5" (49.5 cm)
Back Height:	27.3" (69.3 cm)
Seat Height:	20.3" (51.6 cm)
Overall Width:	33.5" (85.1 cm)
Overall Depth:	72.2" (183.4 cm)
Overall Height (down):	43.7" (111 cm)
Overall Height (raised):	57.7" (146.6 cm)
Max. Vertical Lift:	n/a
Box 1	80CM x 86CM x 83CM
Box 2	n/a
Net weight	53 kg
Gross Weight	59 kg



LIFT CHAIRS (2 MOTORS)

VENUS CH4009M (BROWN) CH4010M (BLUE)

The Venus lift chair is an infinite position lift chair available with brown or blue upholstery. While it comes fitted with rubber pads to prevent floor damage and provide stable positioning, the 4 supplied locking casters can be installed to ease repositioning the chair.

With a Medium/Average size rating, this chair will accommodate people of medium width and average height.

- Infinite position recline

- Four locking swivel casters included

MODEL NAMES	VENUS
Model numbers:	CH4009M (Brown) CH4010M (Blue)
Weight Capacity:	352 lbs (160 kg)
Seat Width:	19.5" (49.5 cm)
Seat Depth:	19.5" (49.5 cm)
Back Height:	29.3" (74.4 cm)
Seat Height:	19.5" (49.5 cm)
Overall Width:	32.4" (82.3 cm)
Overall Depth:	72.9" (185.2 cm)
Overall Height (down):	45.2" (114.8 cm)
Overall Height (raised):	58.5" (148.6 cm)
Max. Vertical Lift:	n/a
Box 1	84CM x 87CM x 68CM
Box 2	85CM x 70CM x 36.5CM
Net weight	47 kg (box1) 6.7 kg (box 2)
Gross Weight	50.3 kg (box1) 8.7 kg (box 2)



JUPITER CH4002M (BLUE) CH4003M (BROWN)

The Jupiter lift chair is an infinite position lift chair available with brown or blue upholstery. While it comes fitted with rubber pads to prevent floor damage and provide stable positioning, the 4 supplied locking casters can be installed to ease repositioning the chair.

With a Medium/Average size rating, this chair will accommodate people of medium width and average height.

- Infinite position recline

- Four locking swivel casters included

MODEL NAMES	JUPITER
Model numbers:	CH4002M (Blue) CH4003M (Brown)
Weight Capacity:	352 lbs (160 kg)
Seat Width:	19.5" (49.5 cm)
Seat Depth:	20.3" (51.6 cm)
Back Height:	28.1" (71.4 cm)
Seat Height:	19.5" (49.5 cm)
Overall Width:	35.1" (89.2 cm)
Overall Depth:	73.3" (186.2 cm)
Overall Height (down):	43.7" (111 cm)
Overall Height (raised):	58.9" (149.6 cm)
Max. Vertical Lift:	n/a
Box 1	85CM x 94CM x 69CM
Box 2	79CM x 83CM x 35CM
Net weight	50 kg (box1) 8.5 kg (box 2)
Gross Weight	54 kg (box1) 10.5 kg (box 2)



LIFT CHAIRS (1 MOTOR)

JUPITER CH4014M (BLUE)

The Single Motor Jupiter lift chair is an three position lift chair available with brown or blue upholstery. While it comes fitted with rubber pads to prevent floor damage and provide stable positioning, the 4 supplied locking casters can be installed to ease repositioning the chair.

Its simple two button pendant control operates the foot and back together with the push of a button. One button for up and one button for down and recline.

With a medium/average size rating, this chair will accommodate people of medium width and average height.

- Three position recline

- Four locking swivel casters included

MODEL NAMES	JUPITER (ONE MOTOR)
Model numbers:	CH4014M (Blue)
Weight Capacity:	352 lbs (160 kg)
Seat Width:	19.5" (49.5 cm)
Seat Depth:	20.3" (51.6 cm)
Back Height:	28.1" (71.4 cm)
Seat Height:	19.5" (49.5 cm)
Overall Width:	35.1" (89.2 cm)
Overall Depth:	73.3" (186.2 cm)
Overall Height (down):	43.7" (111 cm)
Overall Height (raised):	58.9" (149.6 cm)
Max. Vertical Lift:	n/a
Box 1	85CM x 94CM x 69CM
Box 2	79CM x 83CM x 35CM
Net weight	50 kg (box1) 8.5 kg (box 2)
Gross Weight	54 kg (box1) 10.5 kg (box 2)



JUPITER CH4015M (BROWN)

The Single Motor Jupiter lift chair is an three position lift chair available with brown or blue upholstery. While it comes fitted with rubber pads to prevent floor damage and provide stable positioning, the 4 supplied locking casters can be installed to ease repositioning the chair.

Its simple two button pendant control operates the foot and back together with the push of a button. One button for up and one button for down and recline.

With a medium/average size rating, this chair will accommodate people of medium width and average height.

- Three position recline

- Four locking swivel casters included

MODEL NAMES	JUPITER (ONE MOTOR)
Model numbers:	CH4015M (Brown)
Weight Capacity:	352 lbs (160 kg)
Seat Width:	19.5" (49.5 cm)
Seat Depth:	20.3" (51.6 cm)
Back Height:	28.1" (71.4 cm)
Seat Height:	19.5" (49.5 cm)
Overall Width:	35.1" (89.2 cm)
Overall Depth:	73.3" (186.2 cm)
Overall Height (down):	43.7" (111 cm)
Overall Height (raised):	58.9" (149.6 cm)
Max. Vertical Lift:	n/a
Box 1	85CM x 94CM x 69CM
Box 2	79CM x 83CM x 35CM
Net weight	50 kg (box1) 8.5 kg (box 2)
Gross Weight	54 kg (box1) 10.5 kg (box 2)



Planet Series Lift Chairs

Available for Immediate Delivery!


(While Quantities Last)

Models: Earth
Mercury Saturn
Pluto Jupiter
Venus Mars

Features*

Heat & Massage Battery Back Up
Infinite Positioning Lounger Style
Vertical Lift Casters
Incontinence Proof Upholstery

**Features available dependent on model chosen.*



MODEL NAMES	SEAT WIDTH	Seat Depth	Seat Height*	Back Height	Notes
Mercury	19½"	19.1"	19.1"	26.9"	Heat & Massage, Black or Tan Leather Upholstery, Tray, Cup Holder
Venus	19½"	19½"	19½"	29¼"	Brown or Blue
Jupiter	19½"	20¼"	19½"	28.1"	Brown or Blue
Pluto	19½"	19½"	20.3"	27.3"	Brown or Blue
Saturn	21.1"	19½"	19½"	28½"	Brown or Blue, Tray, Cup Holder
Mars	21"	21"	19.9"	26.9"	Vertical Lift, Incontinent Proof Upholstery, Removable Arms
Earth	20¾"	19.1"	19.1"	26.9"	Vertical Lift, Incontinent Proof Upholstery, Removable Arms

EZee Life™



34 Futurity Gate, Unit # 15 Concord Ontario L4K 1S6

GTA: 416-739-8393 | Toll Free: 1-855-235-2400

email: info@factorydirectmedical.com | www.factorydirectmedical.com