

EZee Fold Electric Wheelchairs

	Generation 1	Generation 1	Generation 1	Generation 1	Generation 2	Generation 2	Generation 2	Generation 3	Generation 3	Generation 4	Generation 5
<b>Model Number</b>	<b>CH4050</b>	<b>CH4055</b>	<b>CH4060</b>	<b>CH4064/5</b>	<b>CH4051/4056</b>	<b>CH4053/4054</b>	<b>CH4058/4059</b>	<b>CH4080/ Platinum</b>	<b>CH4085/ Platinum</b>	<b>CH4075</b>	<b>CH4070</b>
<b>Overview</b>	Economy 8" wheels	Economy 8" wheels Removable footrest	Economy 12" Rear wheels	Economy 12" Rear wheels for tall Clients (footboard or swingaway footrests)	Economy 8" wheels/ single battery in side frame (2nd optional) / removal	Two batteries mounted in side frames, LED battery ind., HD forks, flip down, or up, arms	Pro 12" wheels / high weight capacity / 2 underseat batteries	Upgraded suspension, bifold back/leg guards/ inc 2 batteries	Upgraded suspension, bifold back/leg guards/inc 2 batteries	Bariatric	Ultra-light weight
Seat Width	40 cm (15.75")	40 cm (15.75")	45.7 cm (18")	43 cm (17")	40.6 cm (16")	40 cm (15.75")	47 cm (18.5")	42 cm (16.5")	45.7 cm (18")	50.8 cm (20")	46 cm (18.1")
Seat Depth	43 cm (17")	43 cm (17")	40 cm (15.75")	43 cm (17")	40.6 cm (16")	40 cm (15.75")	40 cm (15.75")	42 cm (16.5")	42 cm (16.5")	50.8 cm (20")	40 cm (15.75")
Seat Height	46 cm (18.1")	46 cm (18.1")	44.2 cm (16.6")	50.7 cm (20")	49.5 cm (19.5")	47 cm (18.5")	47 cm (18.5")	51 cm (20")	55 cm (21.7")	54 cm (21.3")	52.5 cm (20.7")
Back Height	45 cm (17.5")	45 cm (17.5")	Adj to 65 cm (25.6")	44.1 cm (17.4")	45.7 cm (18")	43 cm (17")	45.7 cm (18")	45 cm (17.7")	50 cm (19.7")	45.7 cm (18")	39 cm (15.4")
Frame Colour	Silver	Silver	Silver	Silver	Black or Red	Silver or Black	Silver or Orange	Silver	Silver	Silver	Silver
Footrests	Footboard	Swing Away	Footboard	FB or Swing Away	Footboard	Footboard or ELR*	Footboard	Footboard or ELR*	Footboard or ELR*	Footboard	Footboard
Weight Capacity	120 kg (266 lb)	120 kg (266 lb)	160 kg (350 lb)	160 kg (350 lb)	120 kg (266 lb)	120 kg (266 lb)	150 kg (330 lb)	150 kg (330 lb)	150 kg (330 lb)	227 kg (500 lb)	100 kg (220 lb)
Weight of Chair	23.1 kg (51.4 lb)	21.5 kg (47.8 lb)	25.3 kg (56 lb)	25 kgs (55 lb)	23 kg (50.7 lb)	23 kg (50.7 lb)	23 kg (50.7 lb)	27.5 kg (60.5 lb)	33 kg (72.6 lb)	38.2 kg (84 lb)	17 kg (37.5 lb)
Maximum Speed	6.5 kph (4 mph)	6.5 kph (4 mph)	7.2 kph (4.5 mph)	7.2 kph (4.5 mph)	6 kph (3.7 mph)	6 kph (3.7 mph)	6 kph (3.7 mph)	6 kph (3.7 mph)	6 kph (3.7 mph)	8 kph (5 mph)	4.5 KPH (2.8 mph)
Maximum Range	15 km/battery	15 km/battery	15 km/battery	15 km/battery	16 km/battery	20 km/battery	22 km (13.7 miles)	18 km (11.3 miles)	21 km (13.1 miles)	20 km (12.5 miles)	15 km (9.3 Miles)
Maximum Incline	10 degrees	10 degrees	12 degrees	12 degrees	12 degrees	12 degrees	12 degrees	12 degrees	12 degrees	10 degrees	Static Stability: >6°
Front Casters	8"	8"	8"	8"	6"	7"	8"	7"	8"	8"	5"
Drive Wheels	10"	10"	12"	12"	8"	10"	12"	8"	12"	12"	8"
Battery Type	24V Li-ion	24V Li-ion	24V 10A Li-ion	24V 10A Li-ion	24V 10A Li-ion	24V 6A Li-ion	24V 7A Li-ion	12V 16A Li-ion (2)	12V 20A Li-ion (2)	24 V, 20 Ah	24V 10.4A Li-ion
Charging Time	5 Hrs/Battery	5 Hrs/Battery	5 Hrs/Battery	5 Hrs/Battery	8 to 10 hours	6 to 8 hours	8 to 10 hours	8 to 10 hours	8 to 10 hours	8 to 10 hours	5 to 7 hours
											

 EZee Life™

 EZee Fold