

Be Free



WITH **ATTO** SCOOTERS,
NOTHING COMES BETWEEN YOU
AND YOUR DREAMS

movinglife[®]
life changing solutions

www.movinglife.com

Designed to support everyday mobility

Advanced features and functionality



DRIVE MODE
Compact to full size in 10 seconds



SPLIT MODE
Separates into 2 lightweight parts



COMPACT MODE
Folds into a small suitcase size



TROLLEY MODE
Rolls like a suitcase



See the world *and be seen!*



Fits any car trunk,
train and tram



Approved for
onboard flight



Unmatched
travel experience

ATTO SPORT unique features: Breathtaking design, cutting-edge technology

Available on ATTO SPORT, ATTO SPORT MAX, SHABBATTO, SHABBATTO MAX



Front disc brake
for extra safety



Illuminated LCD display /
USB phone charger



Extra wide rear wheel space
offers better stability



Powerful LED headlights
/ Rear hazard light



Shock absorbing
flat free tires



PU leather
seat and grips



FEATURE	ATTO	ATTO SPORT	ATTO SPORT MAX
Max Speed	6.4 kp/h / 4 mph	10 kp/h / 6.2 mph	8 kp/h / 5 mph
Max User Weight	120 kg / 265 lb	120 kg / 265 lb	136 kg / 300 lb
Battery Size	249.6 wh	249.6 wh	249.6 wh
Driving Distance Per Charge	20 km / 12.5 miles	20 km / 12.5 miles	18 km / 11 miles
Weight	Front Half: 26.5 lb / 12 kg Rear Half: 37.5 lb / 17 kg Assembled: 64 lb / 29 kg Battery: 4 lb / 1.8 kg	Front Half: 27.5 lb / 12.5 kg Rear Half: 38.5 lb / 17.5 kg Assembled: 66 lb / 30 kg Battery: 4 lb / 1.8 kg	Front Half: 27.5 lb / 12.5 kg Rear Half: 38.5 lb / 17.5 kg Assembled: 66 lb / 30 kg Battery: 4 lb / 1.8 kg
Cruise Control (Throttle Not Needed)	No	Yes	Yes
Automatic Slowdown	Slopes	Slopes / Turns	Slopes / Turns
Front Disk Brake for Extra Safety	No	Yes	Yes
Display	LED	LCD / Programable: beeps, units, pairing	LCD / Programable: beeps, units, pairing
Tires	Airless PU, flat free tires	Airless PU, shock absorbing, flat free tires	Airless PU, heavy duty, dense, shock-absorbing, flat-free tires
Antitheft Wireless Locking FOB Key	No	Yes	Yes
Embedded Front Headlights	No	Yes	Yes
Color Set	Classic White	Dark Grey	Dark Grey
Seat and Backrest / Grips	PU Seat / Grips	PU Leather Seat / Grips	PU Leather Seat / Grips
Rear Light Lifting Handle	No	Yes	Yes
Eco Mode	No	Limit speed / Increase range	Limit speed / Increase range
Left-hand Option	Optional	Included	Included
Shabbat Remote	No	No	No

SHABBATTO

The World's Best Shabbat Scooter

Based on the ATTO SPORT scooter

Movinglifes' proprietary "Shabbat Control" eliminates the need for expensive external installation

Approved for full operation on Shabbat by the Zomet Institute



FEATURE	SHABBATTO	SHABBATTO MAX
Max Speed	10 kp/h / 6.2 mph	8 kp/h / 5 mph
Max User Weight	120 kg / 265 lb	136 kg / 300 lb
Battery Size	249.6 wh	249.6 wh
Driving Distance Per Charge	Regular use: 20 km / 12.5 miles During Shabbat: 13 km / 8 miles	Regular use: 18 km / 11 miles During Shabbat: 11.5 km / 7 miles
Weight	Front Half: 27.5 lb / 12.5 kg Rear Half: 38.5 lb / 17.5 kg Assembled: 66 lb / 30 kg Battery: 4 lb / 1.8 kg	Front Half: 27.5 lb / 12.5 kg Rear Half: 38.5 lb / 17.5 kg Assembled: 66 lb / 30 kg Battery: 4 lb / 1.8 kg
Cruise Control (Throttle Not Needed)	Yes	Yes
Automatic Slowdown	Slopes / Turns	Slopes / Turns
Front Disk Brake for Extra Safety	Yes	Yes
Display	LCD / Programmable: beeps, units, pairing	LCD / Programmable: beeps, units, pairing
Tires	Airless PU, shock absorbing, flat free tires	Airless PU, heavy duty, dense, shock-absorbing, flat-free tires
Antitheft Wireless Locking FOB Key	Yes	Yes
Embedded Front Headlights	Yes	Yes
Color Set	Dark Grey	Dark Grey
Seat and Backrest / Grips	PU Leather Seat / Grips	PU Leather Seat / Grips
Rear Light Lifting Handle	Yes	Yes
Eco Mode	Limit speed / Increase range	Limit speed / Increase range
Left-hand Option	Included	Included
Shabbat Remote	Yes	Yes

ATTO ACCESSORIES



Flight Kit



*Cane / Crutches
Holders*



*Mobile Phone
Holder*



ATTO Armrests



*ATTO SPORT /
SHABBATTO Armrests*



*Essentials
Pouch*



*Rear Light Lifting
Handle*



*ATTO Left Hand
Throttle Kit*



*Standard
Flight Battery*



Backpack



Seat Cushion



*Carryall
& Cushion*



Cable Kit



Vira-Shield



movinglife[®]

life changing solutions

Movinglife develops life-changing solutions for people with mobility needs. The ATTO is a contemporary line of mobility scooters, that combine smart design with rugged functionality - without compromising on style. ATTO scooters are optimized for travel by any mode of transportation, from compact cars, to trains, cruise ships and airplanes. With its lightweight usability and easy transition between regular, compact and trolley modes, ATTO is the ideal companion to take with you on any journey.

**The world doesn't stop at your corner shop
and neither should you.**