

Model: CH4020 | Black & Grey

Folding effortlessly in seconds, the EZee Travel is designed to fit easily just about anywhere: a small closet, tight elevators, narrow storage places and, of course, the trunk of your car. Take the EZee Travel with you on car or airplane trips and be up and running as soon as you get there.

Folded and stored for the night, it will be ready when you are the next morning. Slide the EZee Travel into the back of a van with virtually no lifting effort and unfold it again in no time. There are no connectors, latches or mechanisms to open when folding or unfolding. It's that easy!

FEATURES:

- Super light weight
- Safe and stable 4 wheel ride
- Ergonomically designed for comfort and functionality
- Folds easily and quickly
- Store in tight places quickly and easily
- Detach the battery pack and carry inside for recharging
- Battery airline friendly for travel
- Compact design
- Rear wheel drive
- Dual articulating front wheels offer added stability
- Automatic electronic braking system



CH4020	
Product Weight (no battery):	61 lb (27.7 kg)
Carrying Capacity:	250 lbs (115 kg)
Shipping Weight:	110 lbs (50 kg)
Product Dimensions:	(L x W x H) 42" x 22" x 36" (106.7 x 56 x 91.4 cm)
Shipping Dimensions:	43 x 23 x 21 in (110 x 59.5 x 52.5 cm)
Folded (vertical):	40 x 22 x 18 in (101 x 57 x 46 cm)
Seat Width:	18 in (46 cm)
Seat Depth:	18.5 in (46 cm)
Back Rest Height:	18 in (46 cm)
Seat Height from floor:	20 in (50.8 cm)
Turning Radius:	29 in (73.7 cm)
Ground Clearance:	3.65 in (9.3 cm)
Range*:	8 miles (13 km)
Battery Weight:	19 lb (8.6 kg)



Range & speed will vary depending on user weight, battery condition and terrain