

Model: CH4035 - Yellow & Black

The EZee Auto-Fold Scooter is an easy to use power folding portable mobility scooter operated with a remote control or manually. The EZee Auto-Fold has a 300 lb. weight capacity and its standard lightweight Lithium battery provides up to 21.7 km on a single charge (range varies depending on user weight, terrain and battery condition).

The EZee Auto-Fold easily fits into the back of even the smallest hatchback. It is also airline friendly and fully compliant and within regulations for Lithium Ion dry cell batteries.

FEATURES:

- Remote Control
- Automatic Powered Folding & Unfolding
- Very Portable
- Excellent for Aircraft & Cruise Ships
- Lithium Battery



CH4035	
Weight Capacity:	136 kg (300 lbs.)
Product Weight:	24 kg (53 lbs.)
Seat Width:	43.2 cm (17")
Seat-to-Floor Height:	35.5 cm (22")
Overall Width:	46 cm (18.1")
Overall Length:	93.7 cm (36.9")
Frame Colour:	Black - Yellow
Frame Type:	Collapsible
Front Wheel Size:	15.2 cm (6")
Rear Wheel Size:	17.8 cm (7")
Wheel Type:	4-Wheel
Turning Radius:	139.7 cm (55")
Max Climbing Angle:	15 degrees
Max Speed:	6 kph (3.7 mph)
Driving Range:	21.7 km (13.5 Miles)
Warranty on Frame:	Three Year Limited
Warranty on Electronics:	One Year Limited



Range & speed varies depending on user weight, terrain and battery condition.