

MERCURY LIFT CHAIR

Mercury CH4011 (Dark Brown) CH4013 (Sand) CH4021 (Black)

The Mercury lift chair is an infinite position lift chair available with dark brown, black or sand coloured genuine top grain leather upholstery. While it comes fitted with rubber pads to prevent floor damage and provide stable positioning, the 4 supplied locking casters can be installed to ease repositioning the chair. Additionally, the Mercury comes equipped with an adjustable tray, a cup holder for user convenience and heat and vibrate.

With a Medium/Short size rating, this chair will accommodate people of average widths and shorter heights.

- Infinite position recline
- Four locking swivel casters included
- Top grain Leather upholstery
- Adjustable tray and cup holder
- Heat and Vibrate



The spring assisted footrest on all models will flip up if it comes in contact with an obstruction when lowering to minimize injury or damage.



Supplied locking swivel casters can be installed on all models to enable easy movement of the chair. Ideal for users who may want to move their chair occasionally and/or rental units.



Our Mercury models include a cup holder, adjustable tray for the convenience of the user.

FEATURES:

- Lounger style that fully supports the legs (no gap between footrest & seat)
- Removable back for easy handling
- Easy to use hand control
- UL / CUL Approved
- Maximum seat elevation with superb stability in full lift position
- Wide construction provides greater side to side stability
- Top steel base for durability & cross supports
- Engineered bushings & washers at pivot points provide quiet smooth operation & durability
- Direct drive activation synchronizes right & left sides for easy opening
- Low voltage motor, 24 volt
- Battery Back Up

MODEL NAMES	MERCURY
Model numbers:	CH4011 (Dark Brown) CH4013(Sand) CH4021(Black)
Weight Capacity:	352 lbs (160 kg)
Seat Width:	19.5" (49.5 cm)
Seat Depth:	19.1" (48.5 cm)
Back Height:	26.9" (68.3 cm)
Seat Height:	19.1" (48.5 cm)
Overall Width:	28.9" (73.4 cm)
Overall Depth:	69.8" (177.3 cm)
Overall Height (down):	42.5" (108 cm)
Overall Height (raised):	55" (139.7 cm)
Max. Vertical Lift:	n/a
Box 1	81cm x 77cm x 65cm
Box 2	76cm x 67cm x 33cm
Net weight	45 kg (box1) 6 kg (box 2)
Gross Weight	47 kg (box1) 8 kg (box 2)

