

DET-306 INFRARED THERMOMETER

Safe & Accurate

The DET-306 infrared forehead thermometer is specifically designed for use on the forehead. It is a no-contact device capable of measuring people's body temperature by detecting the intensity of infrared light emitted from the forehead. It converts the measured heat into a temperature reading displayed on the LCD.

Features:

- Measures forehead temperature
- Contactless operation
- No probe covers needed
- 10 Reading memory
- Beeps when ready
- Large Display
- Dual Scale (°C or °F)
- 1 Second read
- Replaceable Battery



Specifications:

- Temperature Range: Infrared: 32° to 43°C (89.6° to 109.4°F)
- Power: (2) AAA DC3V batteries (batteries NOT included)
- Battery Life: 6000 readings (approximately 1 year)
- Weight (with batteries): 84 grams (3 oz)
- Warranty: 1 year

