EZee Life™

Mobility • Independence • A Better Tomorrow

CARBON FIBRE
ONLY 14.78 kg (32.6 lbs)

LIGHTWEIGHT FOLDABLE ELECTRIC WHEELCHAIR

CH4078



Compact. Lightweight. Travel Ready.

The CH4078 combines a strong carbon fibre frame with a lightweight, foldable design for easy transport and everyday independence.

- Lightweight Carbon Fibre Frame
- Foldable Design for Easy Transport
- Intuitive Joystick Control
- Long-Lasting Lithium Battery
- Electromagnetic Braking System
- Anti-Tip Wheels
- Ideal for Indoor & Flat Surfaces
- Perfect for Travel and Daily Use



User-Friendly Joystick Control

Clear display with speed, battery level and simple controls.



Compact Foldable Design

Folds quickly and easily for storage or transport.

Specifications

Model	CH4078
Frame Material	Carbon Fibre (Lightweight & Strong)
Max User Weight	136 kg (300 lbs)
Net Weight (with battery)	14.78 kg (32.6 lbs)
Gross Weight (boxed)	18 kg (39.7 lbs)
Open Dimensions (L x W x H)	890 x 620 x 930 mm
Folded Dimensions (L x W x H)	650 x 395 x 795 mm
Seat Width	Approx. 450 mm
Drive System	Dual Motors (Left & Right)
Battery Type	24V Lithium Battery (removable)
Charging Input	AC100–240V, 50/60Hz
Maximum Speed	Up to 6 km/h (user adjustable)
Range	Up to 15–20 km*
Braking System	Electromagnetic Braking
Front Wheels	7" (approx.)
Rear Wheels	10" (approx.)

**Range may vary based on user weight, terrain and driving conditions.*

EZee Life™

ADDRESS
HPU Medical Wholesale
15-34 Futurity Gate
Concord, ON L4K 1S6

CALL US
416-739-1267
info@hpumw.com
www.hpumw.com